



FREE WATER TEST REQUEST

We are pleased to offer FREE water testing to determine if pre-treatment is required to ensure proper operation of your UV sterilizer.

INSTRUCTIONS:

Take your water sample after any existing whole house / point-of-entry water treatment equipment you already have installed (if applicable). We require a minimum of 250ml (8 fluid ounces) of water. Use a regular bottled water bottle. Do NOT use a bottle that has previously contained other beverages or liquids, especially juice or soda pop. Fill the bottle completely so there is no air remaining in the container and cap securely. Keep the sample in a dark place until it is mailed. Send the sample along with a completed copy of this form to the address below. If sending from the USA, please use the U.S. Postal Service and describe the contents as "water sample" on the customs declaration. Record the customs value as \$1.

HomePlus Products Inc. ATTN: Water Testing Lab

5-1490 Pearson Place
Kamloops, B.C. V1S 1J9 Canada

Date, Time and Location of Sample: _____

Existing Treatment Equipment: _____

Name: _____

Phone: _____

Email Address: _____

Lab Use Only Below

Test Parameter	Result	Unit	Acceptable Range	Test Accuracy
Turbidity		NTU	<1	+/- 0.1 NTU
UVT, RAW		%	min 75% (>85% preferred)	+/- 0.5%
UVT, FILTERED		%	min 75% (>85% preferred)	+/- 0.5%
pH		-	7.0 to 8.5	+/- 0.4
Total Hardness		mg/l	<120 mg/l (ppm)	+/- 16%
Total Iron (Fe)		mg/l	<0.3 mg/l (ppm)	+/- 3%
Total Manganese (Mn)		mg/l	<0.05 mg/l (ppm)	+/- 6%
TDS		mg/l	<500	+/- 5%

Tested By/Date:

